
NIF Diagnostics

Presentation to NIF Nuclear Astrophysics Workshop



Robert Kauffman Lawrence Livermore National Laboratory

Aug. 28, 2007



NIF will support three national research agendas



The National Ignition Facility
















UCRL-PRES-230665

Capabilities developed for one mission will be also be available for the other users






An extensive set of capabilities is being developed for the National Ignition Program



The National Ignition Facility





The National Ignition Facility

Optical Diagnostics are used for:

- Energy Balance: Incident and scatteredlight
- Laser Plasma Instabilities: Raman & Brillioun Spectra


Full aperture backscatter diagnostic on 2 quads (Q36B, Q31B)

Near backscatter imager diagnostic on 2 quads

Optical Diagnostics








3ω energy diagnostic3ω power diagnostic

FABS and NBI measures backscattered light on quads Q31B and Q36B





FABS (Full Aperture Backscatter System)

- Time resolved (100 psec)
- Spectrally resolved
 SBS .03 nm

 - SRS 6 nm

NBI (Near Backscatter Imager)

- Time integrated SRS and SBS
- Space and time resolved SRS or SBS

Diagnostic capability for laser-plasma interaction experiments

UCRL-PRES-230665

FABS measures backscattered light through the lens transmitted through the final turning mirror



Full Aperture Backscatter System (FABS)



Scattered light through the lens




Schematic showing scattered light in FABS enclosure


UCRL-PRES-230665




The National Ignition Facility



VISAR is a velocity interferometer for measuring reflections from shock waves and other surfaces











Emitted and reflected light (pre heat, shocks)


Optical Diagnostics

> Diagnostic capability for shocked EOS experiments

















The National Ignition Facility




X-Ray
Diagnostics

Hard x-rays
>20 keV

Soft x-rays
.05 to 12 keV
(Hohlraum temperature)

Dante X-ray flux diagnostic for hohlraum drive


Diagnostic capability for X-ray driven experiments





The National Ignition Facility


(1D, 2D, spectroscopy)



2 SXIs (Static X-ray imagers)

NIF-0405-10643r1L6






implosions, spectroscopy, LPI,...


NIF-0405-10643r1L6

There are Diagnostic Instrument Manipulators for inserting diagnostic systems into the chamber






UCRL-PRES-230665



UCRL-PRES-230665








The National Ignition Facility


ARC (Advanced Radiographic Capability)
4 segmented High Irradience Beams

X-Ray Diagnostics



High-energy high-intensity missions define the laser requirements






The petawatt beam path in the target bay will be easily switched from the long-pulse beam path




The National Ignition Facilit



A bright, hard X-ray source enables new diagnostic capabilities on NIF




X-ray Radiography Allows Imaging of Dense, High Z Targets



- Hydrodynamics
- Equation-of-state
- Material Dynamics

X-ray Diffraction Allows Study of Grain Effects

- Material grains




- Material Dynamics
- Phase change kinetics

X-ray Absorption
Spectroscopy is
Greatly Simplified
Implosion backlighter

64 NIF

Beams





- Opacity
- Material dynamics

UCRL-PRES-230665

We have made excellent progress in developing these sources for NIF experiments











The National Ignition Facility



Neutron and gamma diagnostics will be fielded on ignition shots






Detectors will be outside the shield wall with vulnerable components shielded

UCRL-PRES-230665

~30m

Detectors and sensitive electronics will be located outside the 2m-thick shield wall








Neutron and y X-Ray **Nuclear Diagnostics Diagnostics Emissions** neutron and **Backlighter** Hard x-rays Plasma >20 keV **Emitted and reflected light** Soft x-rays (pre heat, shocks) .05 to 12 keV (Hohlraum temperature) **Optical** Scattered ligh **Diagnostics** Incident 300 Diagnostic X-ray Imaging **Support Systems** 5 to 100 keV (1D, 2D, spectroscopy)

NIF-0405-10643r1L8

We successfully fielded ~ half of all the types of diagnostic systems on NIF






Diagnostics returned data 551 out of 577 attempts over a 22 month interval data return was 98% on the last experimental campaign


Diagnostic schedule for NIC







Diagnostic capabilities schedule for science users











Development of common experimental platforms increases shot availability



Example: The NEL Hydrodynamics platform was jointly developed and used for different physics experiments






Most platforms used for NIC and HED are used also for basic science experiments


Capabilities developed for the National Ignition Campaign will be available for science users








Experimental platforms are developed on supporting facilities and demonstrated on NIF

The National Ignition Facility












We are developing the material dynamics platform on for NIF in FY10-FY11











We are developing hydrodynamics platforms for NIF in FY11-FY12






Hydrodynamic Jets



Pure metal foams are being developed for NIF experiments

Hydrodynamic Instabilities



High Mach # RM

We have ~31 types of diagnostic systems planned for NIC through FY11




Diagnostic System		Types	Status	Location
Optical	Full Aperture Backscatter	2	1 Completed, EMP Upgrade Q4FY08; 2nd Q2 FY09	NIF
	Near Backscatter Imager	1	1 Completed, 2nd Q4FY08, 3rd Q3FY09	NIF
	VISAR	1	Completed, relocate	NIF
	Optical Streak Camera	1	Completed	NIF
	Target Illumination	2	1 Each Completed, remainder Q4FY09	NIF
	Pulse Sync	1	Prototyped, Q4FY08	NIF
	Crystal Alignment	1	Prototyped, Q4FY08	NIF
	Dante X-ray Power w/imager	1	Completed	NIF
	Hard X-ray Spectrometer	1	Completed	NIF
	X-ray Streaked Detector	1	Prototyped, 2 Q4FY06	NIF
V Dov	X-ray Gated Detectors	2	Completed	NIF
X-Ray	Front Ends for Detectors	3	4 Prototypes Completed, remainder FY09–FY11	NIF
	Static X-ray Target	2	1 Completed, Gated Q4FY08	NIF
	High Density Target	1	2 Q4FY10, 2 Q4FY11	TITAN, Vulcan
	Hard X-ray Imager	1	Q3FY10	Omega
Neutrons and Gammas	Neutron Time-of-flight w/yield	1	Q3FY10	Nova, Omega
	Neutron Imaging	1	Q3FY10	Nova, Omega, NTS
	Gamma BangTime	1	Q3FY10	Omega, NTS
	Magnetic Recoil Spectroscopy	1	Q3FY10	Omega, NTS
	Activation	1	Q3FY10	Nova, Omega, Z,
	DT Yield	1	Q3FY10	Omega, NTS
Support	Diagnostic Manipulators	1	2 Completed, EMP Upgrades Q4FY07	NIF
	Diagnostic Alignment	1	Completed, EMP Upgrades Q4FY08	NIF
	Chamber Viewing	1	Q1FY07	NIF
	Cross Timing System	1	Exists, expand in time	NIF
Subtotal		31	21 systems tested on NIF	

We successfully fielded 21 out of 31 types of diagnostic systems on NIF Early Light

System Design Requirements (SDR) are developed through Integrated Experimental Teams with stakeholders, reviewed and put under configuration control

The National Ignition Facility

- Elements of SDR
 - Functional description
 - Performance requirements
 - Calibration requirements
 - Environmental Requirements
 - Layout and Utility Interface Requirements
 - Operational Requirements
 - Reliability/Availability/
 Maintainability requirements
 - Quality requirements
- There are typically ~ 60 specifications for each diagnostic system




We are on schedule to complete 90% of all the SDRs this year

UCRL-PRES-230665 4⁴

Example: Dante Soft X-ray Power Diagnostic SDR



Area	SDR	Section	Requirement	
System Description	3.1	What the System does	The soft x-ray power diagnostics measures the absolute soft x-ray emission from a target as a function of time and provides an absolute time history of the radiation temperature inside the hohlraum. The diagnostic also measures "preheat" x-rays in the 2-12 keV range.	
Performance	4.1	Wavelength/Energy Range	The diagnostic shall measure photon energies from 50 eV to 12 keV. The imaging channel should be able to measure photon energies between 500 eV and 3 keV.	
Requirements	4.10	Alignment and Pointing Accuracy	Pointing accuracy in x and y shall be no worse than 10 mm	
Calibrations	5.3	Energy Range of Calibration	The diagnostic shall be calibrated at .05-5 keV with a relative accuracy of ±5% and absolute accuracy of ±10%	
Operational Requirements 8.1 Setup Time After initial installation, set than 15 minutes		After initial installation, setup time shall take no greater than 15 minutes		



We have a requirements plan for diagnostic operations at differing levels of yield



Diagnostic Type	Laser, Energetics, Symmetry, Shock Timing, Ablator Experiments	Convergent shock timing	Ignition Attempts	Ignition
	< 2 1012 DD yield	< 5 e14 TT yield (.5 kJ)	< 1 e16 DT yield (20 kJ)	< 1 e19 DT yield (20 MJ)
Optical	FABS (2), NBI (3), VISAR			
X-Ray	DANTE, FFLEX, GXD, SXD, SXI (2)	SXD (2)	HEXRI	
Neutron &			nTOF, MRS, Activation,	
Particles		nToF	GBT, Protex, NI	Activation
Other	DIM (2), OPAS (2), CTS	DIM (2), OPAS (2), CTS	DIM (2), OPAS (2), CTS	
# of diagnostic				
systems on a				
given shot	~ 16	~ 7	~ 12	~ 1