


# Experimental Platforms on NIF

Warren Hsing  
Presentation to: User Science Workshop on NIF  
February 14, 2012



# We will commission new experimental platforms that can be applied to science experiments

	FY09	FY10	FY11	FY12	FY13	FY14
Implosions & Applications	Hohlraum implosions ★	Direct Drive Exploding pushers ★ ★	Cryogenic DT Implosions ★	Ignition ★		
Plasma physics	★ GasPipe					
Radiation - Hydro		Planer Halfraum ★		★ Spherical: Hohlraum		
High pressure EOS and strength			Hohlraum ICE EOS ★ ★	Hohlraum ICE Strength		
X-ray sources		Low debris 4 KeV source ★	★	Low debris 13 KeV source		
Opacity					★ Long-pulse Hohlraum	




**An integrated suite of capabilities to perform an experiment is termed the experimental “platform”**

**Targets**

**Diagnostics**

**Hohlraum drive**


**Data Analysis**



# Multiple configurations are sometimes needed to address specific physics issues


## Radiation-hydro platform

### Calorimetry configuration





Tr = 200 eV, 10 nsec

### Radiography configuration




Streaked and gated






Once developed, an experimental platform can be customized and applied across a wide variety of physics experiments


### Formation of the Eagle Nebula Pillars



### Non-linear ablation front hydrodynamics




### Formation of Herbig-Haro jets




These experiments can utilize a modified planar-radiation hydrodynamics platform

Turning a concept into a practical experiment requires a suite of capabilities and integrated knowledge of the facility

### EOS concept



### EOS target




# NIF personnel have been identified to work with the current proposals to implement the experiments

## Target capabilities


- Material synthesis
- Capsules
- Micromachining
- Gas fill
- Characterization
- Cryogenics
- Assembly
- Alignment & target fiducials
- Unconverted light & shields
- Debris
- Materials restrictions

## Foam synthesis & micromachining




Low-Z, high-Z and pure metal foams can be synthesized

## Gas Fill



Gas pipes, hohlraums & capsules

## Characterization




X-ray radiographs

## Cryogenics, Unconverted light shields, Assembly

Shock Timing Configuration

Cryo Thermo-Mechanical Package






# NIF personnel have been identified to work with the current proposals to implement the experiments

## Diagnostic capabilities

- Modifications / new design design, build, qualify
- Lines-of-sight
- Detector configuration
- Filters
- X-ray backlighters
- Debris and damage
- Alignment
- Calibration
- Cross-timing
- Materials compatibility
- Unconverted / scattered light
- Activation / contamination





- X-ray diagnostics
  - Time-resolved, time-integrated detectors, imaging & spectral snouts
  - Soft x-ray flux
- Optical diagnostics (Backscatter, VISAR)
- Nuclear diagnostics (nToF, activation, MRS, NI, WRF, Gamma BT, Radchem)

# NIF personnel have been identified to work with the current proposals to implement the experiments


## Laser drive capabilities

- Hohlraum vs direct drive
- Pulse shape
- Beam pointing and focus
- Beam smoothing
- Wavelength separation
- Drive characterization
- Drive symmetry
- Backscatter & hot electrons
- Unconverted & counter-propagating light
- Optics lifetime


## Beam pointing, focus, pulse shape




## $1\omega$ , $2\omega$ unconverted light



## Drive Temperature



## Beam smoothing



CPPs, polarization,  
SSD bandwidth

# NIF personnel have been identified to work with the current proposals to implement the experiments


## Simulations

- Hohlraum / direct drive designs
- Backscatter assessment
- Preheat assessment
- Material database
- X-ray backlighter designs
- Debris, EMP


## Data Analysis

- Instrument correction
- Data unfold routines
- Film/image plate scanning
- Data archival
- Data access


## Hohlraums



## Backlighter debris



## Instrument correction






## It is beneficial to use existing platforms when applicable

Typical lead times prior to a NIF experiment			
Capability	Identical to existing platform	Small modification to existing platform	New capability
Targets	1 – 6 months	6 months – 1 year	> 1 year
Laser drive	~ 1 month	1 – 6 months	> 6 months
Diagnostics	~ 1 – 3 months	3 – 6 months	> 6 months
Data analysis	Exists	1 – 3 months	> 3 months


**Identification of key scientific experiments can drive new capabilities required**

# The National Ignition Campaign uses 4 main platforms




# NIC platforms can be applied to fundamental science

## Convergent Ablator Platform



**Potential applications:**  
Gbar EOS, convergent hydrodynamics (RT, RM, KH)

## Keyhole Platform




**Potential applications:**  
Convergent shock physics, RM instability



# NIC platforms can be applied to fundamental science


## Symcap Platform




**Potential applications:**  
Hohlraum energetics/Laser  
Plasma Instabilities, drive  
characterization for experiments,

## Exploding Pusher Platform

Neutron time-of-flight  
Gaseous radchem  
Solid radchem  
Yield





spectroscopy



**SiO<sub>2</sub> or CH shell  
+ gas dopants  
+ shell dopants**


**Potential applications:**  
NLTE atomic physics,  
spectroscopy, nuclear physics  
w/ radchem

# Planar X-ray driven platforms can reach 10's of Mbars pressure quasi-isentropically or 100's Mbars Hugoniot



## Existing platforms continually evolve


**Planar Hydrodynamics Platforms  
now achieve >300 eV drive**



**Vacuum Half  
hohlraum  
350 eV**

**Gas filled  
hohlraum  
>300 eV**

**X-ray sources will use Fe foam  
filled tubes and Kr gas fill**




**Gas filled tubes**




# Experiments have demonstrated a high degree of reproducibility in the laser and X-ray drive

Laser energy vs request is typically within ~ 3%



Example of pulse shape changes during a NIC campaign


Reproducibility in X-ray flux was measured to be ~ 5%




Example of measurements taken on 3 shots over a 9 month period

# Diagnostics and targets are calibrated or metrologized to a high level of accuracy

**DANTE has ~ 5000 calibrations taken to provide x-ray flux with an absolute accuracy of ~ 7%**




**X-ray radiographs measure uniformity and density-length of completed target to ~ 3%**



**Together with the laser, allows highly precise measurements to be taken on NIF**

## We enhance shot opportunities on NIF by minimizing changes in facility configurations




### Diagnostic Configuration: 3 Yield w/neutron imaging



- A facility configuration is defined to a set of diagnostics, optics, classification and yield configurations that require multiple shifts to change
- Examples of capabilities that require multiple shifts to reconfigure
  - Exchanging instruments in a DIM
  - Yield Experiment  $> 1e14$
  - $>8$  CPPs exchanged
  - Classified diagnostics
  - Special diagnostic reconfigurations
- Examples of capabilities that can be modified during the shot cycle
  - Imaging front end for a gated imager
  - Attenuators and filters on DANTE
  - Sweep speed and timing for a gated instrument
  - These do not define unique configurations

The integrated NIF plan is based on a minimal set of configurations that both meet user needs and optimize facility availability




	Jan	Feb	Mar	Apr	May	Jun	July	Aug	Sept
Configurations	1	2	3	4	5	3	6	3	5

Diagnostic Configuration 1: Compton Radiography	Diagnostic Configuration 2: Conv ablator 90-315	Diagnostic Configuration 3: DT w/neutron imager
		
Compatible with: Symcap, keyhole, EOS, EP, DIME, Pleiades, Fanbolt, viewfactor, Abl RT, EPEC-A	Compatible with: EP, Fanbolt, Abl RT, XRSA, EPEC-A; Alternate config for Symcap, DIME	Compatible with: Symcap, EP, THD, DIME, Mix symcap, Pleiades, Fanbolt, Abl RT, XRSA, EPEC-A



The integrated NIF plan is based on a minimal set of configurations that both meet user needs and optimize facility availability

	Jan	Feb	Mar	Apr	May	Jun	July	Aug	Sept
Configurations	1	2	3	4	5	3	6	3	5

Diagnostic Configuration 4: HED EOS	Diagnostic Configuration 5: Largest number of platforms	Diagnostic Configuration 6: Conv Abl 90-78
		
Compatible with: Matl Strength Drive	Compatible with: Symcap, Keyhole, EP, DT w/o NI, DIME, Pleiades, Fanbolt, Viewfactor, EOS, Abl RT, XRSA, EPEC-A, Mix Symcap	Compatible with: Symcap, Keyhole, EP, DIME, Pleiades, Fanbolt, Mix sym, Viewfactor, EOS, Gbar EOS, Abl RT, XRSA, EPEC-A

# NIF is operational, conducting experiments and acquiring great data



