X-Ray Diffraction on NIF NIF Users Group Meeting LLNL, 2/15/2012

Jon Eggert (LLNL)

Justin Wark (Univ. Oxford)

This work performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344.




Team

- Lawrence Livermore Laboratory
 - J.H. Eggert , J.R. Rygg, , J.A. Hawreliak, R. Smith, D. Swift, D. Braun, P. Celliers, D.H. Kalantar, G. Collins
- University of Oxford
 - J. Wark, A. Higginbotham, M. Suggit, G. Mogni
- National Research Council of Canada
 - D. Klug, Y. Yao
- Washington State University
 - C.S. Yoo
- AWE Aldermaston
 - N. Park
- Osaka University
 - K. Shigemori, K. Shimizu
- University of Rochester
 - Tom Boehly
- Additional Collaborators / Consultants
 - A. Lazicki, F. Coppari, D. Hicks, Y. Ping, D. Swift
 - Andrew Comley, Brian Maddox, Hye-Sook Park, and Bruce Remington
 - M. McMahon, T. Duffy




The National Ignition Facility (NIF) is currently a 192 beam, 1.6 MJ laser



NIF Ramp-Compression Experiments have already made the relevant exo-planet pressure range from 1 to 50 Mbar accessible.


We measured stress-density up to 5 TPa on NIF.

187 transiting exoplanets as of December 29, 2011, http://exoplanet.eu/




Our NIF experiments have demonstrated that we can access the relevant pressure composition region for exo-planet interiors.

We Proposed to Study Carbon Phases by X-Ray Diffraction on NIF



Just a few years ago, ultra-high pressure phase diagrams for materials were very "simple"






New experiments and theories point out surprising and decidedly complex behavior at the highest pressures considered.


Traditional view: All materials become simple at high pressure appears to be incorrect!

"... what the present results most assuredly demonstrate is the importance of pressure in revealing the limitations of previously hallowed models of solids"

—Neil Ashcroft (2009).




High pressures phases of aluminum are also predicted to be complex



Diffraction on NIF

Recent metadynamics survey of carbon proposed a dynamic pathway among multiple phases



Canales, Pickard, Needs, Phys. Rev. Lett. 108, 045704 (2012)



A new paradigm. Really?


Are we really about to witness a true paradigm shift in extreme compressed-matter physics?

"Only as experiment and tentative theory are together articulated to a match does the discovery emerge and the theory become a paradigm"

-p 61

"Further development ordinarily calls for the construction of elaborate equipment, the development of an esoteric vocabulary and skills, and a refinement of concepts..."

-р 64





Facilitating new Paradigms in Extreme Compression and High Energy Density Science



We need to develop diagnostics and techniques to explore this new regime of highly compressed-matter science.



X-Ray Diffraction:

Understand the phase diagram / EOS / strength / texture of materials to 10's of Mbar

Strategy and physics goals:

Powder diffraction

Begin with diamond

Continue with metals etc.

Explore phase diagrams

Develop liquid diffraction


Reduce background / improve resolution




Powder x-ray diffraction of rolled foils on the Omega laser





We performed high-pressure x-ray diffraction on tantalum at the Omega laser






Diffraction data quality is roughly where DAC diffraction was in the '80s. We need to make similar strides.







We determine stress by backward integration of diamond freesurface velocity



Tantalum diffraction on the Omega laser





Results:

- High-quality data at moderate pressure
- Extension of the bcc equation of state
- Indication of possible phase transition above 300 GPa 600 GPa

Further work:
Solve technical difficulties


associated with diffraction above


Lessons learned on the Omega laser:


- Samples are cool enough to do crystalline diffraction up to near
 1 Tpa, far above Hugoniot melt pressure
- Small number of reflections available is a major limitation in structure determination
- X-ray background is a primary concern
- Above ~300 GPa texturing is ubiquitous (for both single- and poly-crystalline materials)



Texture and the Ewald Sphere Construction












Texture and the Ewald Sphere Construction








Texture and the Ewald Sphere Construction








White light x-ray diffraction of single crystals on the Omega laser

Experiments performed by Andrew Comley, Brian Maddox, Jim Hawreliak, Hye-Sook Park, and Bruce Remington




We probe shocked Ta (100) crystals in-situ using white-light Laue x-ray diffraction at the OMEGA laser facility



Thanks to Andrew Comley, Brian Maddox, Hye-Sook Park, and Bruce Remington



(difference in strains in shocked and transverse directions) of the compressed unit cell



- Von Mises Stress (Strength) = 2C' x strain anisotropy
- C' = $(C_{11} C_{12}) / 2$ where C_{11} and C_{12} are elastic constants

Thanks to Andrew Comley, Brian Maddox, Hye-Sook Park, and Bruce Remington





Guiding principles for fielding x-ray diffraction on NIF:

- Design and field a diagnostic useful to a wide range of experiments
 - Powder diffraction
 - White-light single crystal Laue
 - EXAFS
- Employ successful designs from Omega
- Enable both direct and indirect drive configurations
- Explore advanced concept possibilities
- Concentrate on reducing background, increasing resolution, and increasing the number of reflections observed





Pathway to the NIF

Diagnostic development

- Responsible Scientist: Ray Smith
- Responsible Individual: John Dzenitis
- Qualify diagnostic for
 - Debris
 - Survivability
- Determine
 - Optimum shielding for hohlraum drive
 - Optimum backlighter energy
- CDR planned in 2 months


Shot Plan


- 1. Diagnostic Damage assessment (½ Energy 500 GPa)
- 2. Noise Level measurement of Image plates
- 3. Low (1 1.2 TPa)
- 4. Middle (1.7 2.0 TPa)
- 5. High (2.5 3.0 Tpa)





We will combine the XRD capability from Omega and the successful drive already used on the NIF











The NIF targets will draw on Omega design and experience







Using target components which will not contribute to the diffraction signal we ensure we probe a limited temporal and spatial region at the peak pressure






Diagnostic configuration





Advanced Designs: Soller Slits.



Advanced Designs: Energy / Angle Dispersive Diffraction

Using a broadband x-ray source and a fixed location energy dispersive detector we can resolve different lattice planes at different energies (λ = 2d sin θ where θ is fixed and λ varies as d)

The example geometry shown below is similar to a pinhole camera which images the sample material onto a CCD. By filtering and use of small pinhole apertures the single

level can be reduced to a single photon counting level where the CCD can provide spectral resolution Broadband X Plastic/Al shielding from (Implosion or shor Source to CCD High-Z pinhole camera housing high energy particles d X-rays and louresence Broadband X-ray Source CCD (Implosion or shortpulse on foil) Source to CCD Diffracted X-rays and Sample flouresence 1 mm thick. amera housing might energy particles Bro 1 mm thick. (Implosi Source to CCD **VISAR** 0.5 mm diameter 1 mm thick, High-Z aperature 1 mm diameter Diffracted X-rays and High-Z aperature Sample flouresence aperature flourscence from first aperature.


This design would allow coincident EXAFS measurements








Need in situ powder X-ray diffraction to make direct measurement of atomic structure






Each possible structure has a unique x-ray diffraction signal that can be measured

Powder XRD can distinguish all of the proposed carbon phases.

Extended X-ray Absorption Fine Structure (EXAFS) at the Omega Laser

temperature



Hicks,Ping



Quantitative structural data by single-shot EXAFS on the Omega Laser

